



Welch Allyn Pocket LED Otoscopes, Ophthalmoscopes and Sets

Maximised performance and durability in a compact set that travels easily with you from one patient to the next.



Examine your patients in style


with high-performance diagnostic tools from Welch Allyn



Four different device colours to choose from:



Onyx



Plum



Vanilla



Blueberry

Customise with our accessory kit* to mix and match your look.



Welch Allyn Pocket LED Handle Bumpers and Otoscope Window Bumpers

Also available as sets:

Welch Allyn
Pocket LED Sets
with Soft Case



Welch Allyn
Pocket Plus LED Sets
with Hard Case



**Pocket Plus LED set comes standard with accessories in all four colours.
Pocket Plus LED otoscope comes standard with matching colour and vanilla accessories.
Pocket Plus LED ophthalmoscope comes standard with matching colour and vanilla handle bumpers.*

Welch Allyn Pocket LED Otoscope

- The most durable, pocket-style otoscope on the market today, based on rigorous drop test versus leading competitors
- Portable, to help you meet the challenges you face when treating patients
- Lightest-weight Welch Allyn otoscope, for the comfort and control you need
- Fibre optics provide cool light with no reflections or obstructions



Pocket Plus LED Otoscope

All the same features as the Pocket LED otoscope, plus:

- Six lumens for 2X the light intensity
- 5-year warranty

Welch Allyn Pocket LED Ophthalmoscope



- Featuring Welch Allyn SureColor™ LED lamp technology, designed to provide the exact light you need for fundoscopic examination
- Lightest-weight Welch Allyn ophthalmoscope for portability and ease of use
- One of the most durable, pocket-style ophthalmoscopes on the market; designed to be a reliable tool for your practice




Pocket Plus LED Ophthalmoscope

All the same features as the Pocket LED ophthalmoscope, plus:

- Neutral density light filter switch, to help virtually eliminate corneal glare
- Soft rubber brow rest, for enhanced comfort during an exam
- 5-year warranty

Welch Allyn SureColor LED Technology



The Welch Allyn LED Difference

With the right balance of colour temperature, colour appearance (or CRI) and brightness, we've optimized LEDs for EENT instruments, providing the benefits of LED light without changing the view you are proficient at for diagnosis.

- Welch Allyn LED Lamps are rated for up to 20,000 hours of use
- Welch Allyn LED lamps offer true tissue colour, clarity and definitions of the ear canal and fundus
- No lamp replacement required

To learn more about Pocket LED and Pocket Plus LED Oscopes, Ophthalmoscopes and Sets, contact your Welch Allyn representative, or visit welchallyn.co.uk

huisman medisch

WelchAllyn®